

For Research Use Only

AP4B1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84710-2-PBS



Basic Information

Catalog Number: 84710-2-PBS	GenBank Accession Number: BC014146	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10717	CloneNo.: 242121G11
Source: Rabbit	UNIPROT ID: Q9Y6B7	
Isotype: IgG	Full Name: adaptor-related protein complex 4, beta 1 subunit	
Immunogen Catalog Number: AG9348	Calculated MW: 739 aa, 83 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01515-1, AP4B1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84710-2-PBS. Detection antibody: 84710-1-PBS. Standard: Ag9348. Range: 0.781-100 ng/mL.