

For Research Use Only

WASP Recombinant antibody

Catalog Number:84714-2-RR



Basic Information

Catalog Number:

84714-2-RR

Size:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1444

GenBank Accession Number:

BC012738

GeneID (NCBI):

7454

UNIPROT ID:

P42768

Full Name:

Wiskott-Aldrich syndrome (eczema-thrombocytopenia)

Calculated MW:

53 kDa

Observed MW:

58 kDa

Purification Method:

Protein A purification

CloneNo.:

242074E4

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Jurkat cells, TF-1 cells, HL-60 cells, Raji cells

Background Information

Wiskott-Aldrich Syndrome protein (WASP) regulates actin cytoskeleton in hematopoietic cells and defects in WASP cause Wiskott-Aldrich Syndrome (WAS), an X-linked immunodeficiency and autoimmunity disorder of childhood.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84714-2-RR (WASP antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 84714-2-RR against Human WASP were performed. The affinity constant is 1.88 nM.