

For Research Use Only

FOSL2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84718-4-PBS



Basic Information

Catalog Number: 84718-4-PBS	GenBank Accession Number: BC008899	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2355	CloneNo.: 242225F7
Source: Rabbit	UNIPROT ID: P15408	
Isotype: IgG	Full Name: FOS-like antigen 2	
Immunogen Catalog Number: AG8397	Calculated MW: 326 aa, 35 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01531-1, FOSL2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84718-4-PBS. Detection antibody: 84718-1-PBS. Standard: Ag8397. Range: 0.313-40 ng/mL