

For Research Use Only

NOTCH3 Recombinant antibody, PBS Only (Capture)

Catalog Number:84721-1-PBS



Basic Information

Catalog Number:

84721-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_000435

GeneID (NCBI):

4854

UNIPROT ID:

Q9UM47

Full Name:

Notch homolog 3 (Drosophila)

Calculated MW:

244 kDa

Purification Method:

Protein A purification

CloneNo.:

242267D1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01506-1, NOTCH3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84721-1-PBS. Detection antibody: 84721-2-PBS. Standard: Ag34533. Range: 0.625-80 ng/mL.