

For Research Use Only

ARMC9 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84724-2-PBS



Basic Information

Catalog Number: 84724-2-PBS	GenBank Accession Number: BC004514	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 80210	CloneNo.: 242346D5
Source: Rabbit	UNIPROT ID: Q7Z3E5	
Isotype: IgG	Full Name: armadillo repeat containing 9	
Immunogen Catalog Number: AG27511	Calculated MW: 817 aa, 92 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01533-1, ARMC9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84724-1-PBS. Detection antibody: 84724-2-PBS. Standard: Ag27511. Range: 0.625-80 ng/mL