

For Research Use Only

Mouse MECR Recombinant antibody, PBS Only (Detector)

Catalog Number:84729-1-PBS



Basic Information

Catalog Number: 84729-1-PBS	GenBank Accession Number: BC003864	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 26922	CloneNo.: 242151D9
Source: Rabbit	UNIPROT ID: Q9DCS3	
Isotype: IgG	Full Name: mitochondrial trans-2-enoyl-CoA reductase	
Immunogen Catalog Number: AG0465	Calculated MW: 40 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01508-2, MOUSE MECR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84729-2-PBS. Detection antibody: 84729-1-PBS. Standard: Ag0465. Range: 0.625-40 ng/mL.