

For Research Use Only

Mouse MECR Recombinant antibody, PBS Only (Capture)

Catalog Number:84729-2-PBS



Basic Information

Catalog Number:

84729-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0465

GenBank Accession Number:

BC003864

GeneID (NCBI):

26922

UNIPROT ID:

Q9DCS3

Full Name:

mitochondrial trans-2-enoyl-CoA
reductase

Calculated MW:

40 kDa

Purification Method:

Protein A purification

CloneNo.:

242151D11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01508-1, MOUSE MECR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84729-2-PBS. Detection antibody: 84729-3-PBS. Standard: Ag0465. Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP01508-2, MOUSE MECR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84729-2-PBS. Detection antibody: 84729-1-PBS. Standard: Ag0465. Range: 0.625-40 ng/mL.