

For Research Use Only

PER1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84730-1-PBS



Basic Information

Catalog Number: 84730-1-PBS	GenBank Accession Number: BC028207	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5187	CloneNo.: 242096A4
Source: Rabbit	UNIPROT ID: O15534	
Isotype: IgG	Full Name: period homolog 1 (Drosophila)	
Immunogen Catalog Number: AG4249	Calculated MW: 1290 aa, 136 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01530-1, PER1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84730-1-PBS. Detection antibody: 84730-3-PBS. Standard: Ag4249. Range: 0.313-40 ng/mL