

For Research Use Only

# TXNDC5 Recombinant antibody, PBS Only (Capture)

Catalog Number:84756-4-PBS



## Basic Information

Catalog Number:

84756-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13965

GenBank Accession Number:

BC001199

GeneID (NCBI):

81567

UNIPROT ID:

Q8NBS9

Full Name:

thioredoxin domain containing 5  
(endoplasmic reticulum)

Calculated MW:

432 aa, 48 kDa

Purification Method:

Protein A purification

CloneNo.:

242195D2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01536-1, TXNDC5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84756-4-PBS. Detection antibody: 84756-1-PBS. Standard: Ag13965. Range: 1.56-100 ng/mL.