

For Research Use Only

CACNB3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84762-2-PBS



Basic Information

Catalog Number: 84762-2-PBS	GenBank Accession Number: NM_000725.3	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 784	CloneNo.: 242088D6
Source: Rabbit	UNIPROT ID: P54284	
Isotype: IgG	Full Name: calcium channel, voltage-dependent, beta 3 subunit	
Immunogen Catalog Number: AG35651	Calculated MW: 55kDa, 484aa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01545-1, CACNB3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84762-3-PBS. Detection antibody: 84762-2-PBS. Standard: Ag35651. Range: 0.156-20 ng/mL.