

For Research Use Only

HMGCL Recombinant antibody, PBS Only (Capture)

Catalog Number: 84769-2-PBS



Basic Information

Catalog Number: 84769-2-PBS	GenBank Accession Number: BC010570	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3155	CloneNo.: 242224H1
Source: Rabbit	UNIPROT ID: P35914	
Isotype: IgG	Full Name: 3-hydroxymethyl-3-methylglutaryl-Coenzyme A lyase	
Immunogen Catalog Number: AG10426	Calculated MW: 325 aa, 34 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01555-2, HMGCL Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84769-2-PBS. Detection antibody: 84769-1-PBS. Standard: Ag10426. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01555-1, HMGCL Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84769-2-PBS. Detection antibody: 84769-3-PBS. Standard: Ag10426. Range: 0.313-40 ng/mL