

For Research Use Only

MAGI3 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84771-2-PBS



Basic Information

Catalog Number: 84771-2-PBS	GenBank Accession Number: BC130409	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 260425	CloneNo.: 242212C5
Source: Rabbit	UNIPROT ID: Q5TCQ9	
Isotype: IgG	Full Name: membrane associated guanylate kinase, WW and PDZ domain containing 3	
Immunogen Catalog Number: AG18394	Calculated MW: 1481 aa, 163 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01559-1, MAGI3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84771-2-PBS. Detection antibody: 84771-1-PBS. Standard: Ag18394. Range: 0.313-20 ng/mL



Cytometric bead array standard curve of MP01559-2, MAGI3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84771-3-PBS. Detection antibody: 84771-2-PBS. Standard: Ag18394. Range: 0.313-20 ng/mL