

For Research Use Only

# ACSS2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84774-2-PBS



## Basic Information

Catalog Number:

84774-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9121

GenBank Accession Number:

BC010141

GeneID (NCBI):

55902

UNIPROT ID:

Q9NR19

Full Name:

acyl-CoA synthetase short-chain  
family member 2

Calculated MW:

701 aa, 79 kDa

Purification Method:

Protein A purification

CloneNo.:

242253C12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01551-2, ACS52 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84774-4-PBS. Detection antibody: 84774-2-PBS. Standard: Ag9121. Range: 0.781-100 ng/mL.