

For Research Use Only

RAB21 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84776-3-PBS



Basic Information

Catalog Number: 84776-3-PBS	GenBank Accession Number: NM_014999	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 23011	CloneNo.: 242117C6
Source: Rabbit	UNIPROT ID: Q9UL25	
Isotype: IgG	Full Name: RAB21, member RAS oncogene family	
Immunogen Catalog Number: AG35291	Calculated MW: 24 kDa	

Applications

Tested Applications:
Cytometric bead array, Sandwich ELISA, Indirect ELISA,
Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Sandwich ELISA standard curve of MP01552-1, Human RAB21 Recombinant Matched Antibody Pair - PBS only. 84776-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag35291. 84776-4-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL



The mean RAB21 concentration was determined to be 10.3 ng/mL in A549 cell extract based on a 2.5 mg/mL extract load, 9.5 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load and 5.5 ng/mL in HepG2 cell extract based on a 2.3 mg/mL extract load.



Cytometric bead array standard curve of MP01552-1, RAB21 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84776-3-PBS. Detection antibody: 84776-4-PBS. Standard: Ag35291. Range: 0.625-80 ng/mL