

For Research Use Only

LRRFIP1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84778-2-PBS



Basic Information

Catalog Number:

84778-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34760

GenBank Accession Number:

BC108914

GeneID (NCBI):

9208

UNIPROT ID:

Q32MZ4

Full Name:

leucine rich repeat (in FII)
interacting protein 1

Purification Method:

Protein A purification

CloneNo.:

242161F10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01600-1, LRRFIP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84778-1-PBS. Detection antibody: 84778-2-PBS. Standard: Ag34760. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP01600-2, LRRFIP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84778-3-PBS. Detection antibody: 84778-2-PBS. Standard: Ag34760. Range: 0.313-40 ng/mL.