

For Research Use Only

LRRFIP1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84778-3-PBS



Basic Information

| | | |
|---|---|---|
| Catalog Number: 84778-3-PBS | GenBank Accession Number: BC108914 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 9208 | CloneNo.: 242161F4 |
| Source: Rabbit | UNIPROT ID: Q32MZ4 | |
| Isotype: IgG | Full Name: leucine rich repeat (in FLII) interacting protein 1 | |
| Immunogen Catalog Number: AG34760 | | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01600-2, LRRFIP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84778-3-PBS. Detection antibody: 84778-2-PBS. Standard: Ag34760. Range: 0.313-40 ng/mL.