

For Research Use Only

BIRC3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84779-3-PBS



Basic Information

Catalog Number: 84779-3-PBS	GenBank Accession Number: BC027485	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 330	CloneNo.: 242214G1
Source: Rabbit	UNIPROT ID: Q13489	
Isotype: IgG	Full Name: baculoviral IAP repeat-containing 3	
Immunogen Catalog Number: AG20500	Calculated MW: 604 aa, 68 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01556-1, BIRC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84779-3-PBS. Detection antibody: 84779-1-PBS. Standard: Ag20500. Range: 0.313-20 ng/mL



Cytometric bead array standard curve of MP01556-2, BIRC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84779-3-PBS. Detection antibody: 84779-2-PBS. Standard: Ag20500. Range: 0.313-20 ng/mL