

For Research Use Only

FGD1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84783-1-PBS



Basic Information

Catalog Number: 84783-1-PBS	GenBank Accession Number: BC034530	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2245	CloneNo.: 242314C2
Source: Rabbit	UNIPROT ID: P98174	
Isotype: IgG	Full Name: FYVE, RhoGEF and PH domain containing 1	
Immunogen Catalog Number: AG3861	Calculated MW: 901 aa, 91 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01538-1, FGD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84783-2-PBS. Detection antibody: 84783-1-PBS. Standard: Ag3861. Range: 6.25-200 ng/mL