

For Research Use Only

EMR1/ADGRE1/(F4/80) Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84786-1-PBS



Basic Information

Catalog Number:

84786-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001974.5

GeneID (NCBI):

2015

UNIPROT ID:

Q14246-1

Full Name:

egf-like module containing, mucin-like, hormone receptor-like 1

Calculated MW:

98kDa

Purification Method:

Protein A purification

CloneNo.:

241999B5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01597-1, EMR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84786-3-PBS. Detection antibody: 84786-1-PBS. Standard: Eg1716. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP01597-2, EMR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84786-2-PBS. Detection antibody: 84786-1-PBS. Standard: Eg1716. Range: 0.78-100 ng/mL