

For Research Use Only

MCT4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84797-2-PBS



Basic Information

Catalog Number: 84797-2-PBS	GenBank Accession Number: BC112269	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9123	CloneNo.: 242296C3
Source: Rabbit	UNIPROT ID: O15427	
Isotype: IgG	Full Name: solute carrier family 16, member 3 (monocarboxylic acid transporter 4)	
Immunogen Catalog Number: AG18788	Calculated MW: 465 aa, 49 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01590-2, SLC16A3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84797-2-PBS. Detection antibody: 84797-3-PBS. Standard: Ag18788. Range: 0.781-50 ng/mL.



Cytometric bead array standard curve of MP01590-1, SLC16A3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84797-2-PBS. Detection antibody: 84797-1-PBS. Standard: Ag18788. Range: 0.781-50 ng/mL.