

For Research Use Only

MCT4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84797-3-PBS



Basic Information

Catalog Number:

84797-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18788

GenBank Accession Number:

BC112269

GeneID (NCBI):

9123

UNIPROT ID:

O15427

Full Name:

solute carrier family 16, member 3
(monocarboxylic acid transporter 4)

Calculated MW:

465 aa, 49 kDa

Purification Method:

Protein A purification

CloneNo.:

242296G7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01590-2, SLC16A3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84797-2-PBS. Detection antibody: 84797-3-PBS. Standard: Ag18788. Range: 0.781-50 ng/mL.