

For Research Use Only

# ATIC Recombinant antibody, PBS Only (Capture)

Catalog Number: 84798-4-PBS



## Basic Information

Catalog Number:

84798-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1052

GenBank Accession Number:

BC008879

GeneID (NCBI):

471

UNIPROT ID:

P31939

Full Name:

5-aminoimidazole-4-carboxamide  
ribonucleotide formyltransferase/IMP  
cyclohydrolase

Calculated MW:

65 kDa

Purification Method:

Protein A purification

CloneNo.:

242324G4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01575-2, ATIC Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84798-4-PBS. Detection antibody: 84798-2-PBS. Standard: Ag1052. Range: 1.563-100 ng/mL.