

For Research Use Only

Mouse CD9 Recombinant antibody, PBS Only (Detector)

Catalog Number:84801-1-PBS



Basic Information

Catalog Number:

84801-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_007657.3

GeneID (NCBI):

12527

UNIPROT ID:

P40240

Full Name:

CD9 antigen

Calculated MW:

25kDa

Purification Method:

Protein A purification

CloneNo.:

242432B7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01578-2, MOUSE CD9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84801-2-PBS. Detection antibody: 84801-1-PBS. Standard: Eg1397. Range: 0.313-40 ng/mL.