

For Research Use Only

EGF Recombinant antibody, PBS Only (Detector)

Catalog Number: 84802-1-PBS



Basic Information

Catalog Number: 84802-1-PBS	GenBank Accession Number: BC093731	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1950	CloneNo.: 242278A3
Source: Rabbit	UNIPROT ID: P01133	
Isotype: IgG	Full Name: epidermal growth factor (beta-urogastrone)	
Immunogen Catalog Number: HZ-1326	Calculated MW: 1207 aa, 134 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01595-1, EGF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84802-2-PBS. Detection antibody: 84802-1-PBS. Standard: HZ-1326. Range: 3.125-200 ng/mL.