

For Research Use Only

ZDHHC5 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84803-1-PBS



Basic Information

Catalog Number: 84803-1-PBS	GenBank Accession Number: BC026967	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 25921	CloneNo.: 242192B8
Source: Rabbit	UNIPROT ID: Q9C0B5	
Isotype: IgG	Full Name: zinc finger, DHHC-type containing 5	
Immunogen Catalog Number: AG15022	Calculated MW: 715 aa, 78 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01576-2, ZDHHC5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84803-1-PBS. Detection antibody: 84803-2-PBS. Standard: Ag15022. Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP01576-1, ZDHHC5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84803-1-PBS. Detection antibody: 84803-3-PBS. Standard: Ag15022. Range: 0.625-80 ng/mL.