

For Research Use Only

# ZDHHC5 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84803-2-PBS



## Basic Information

|                                             |                                                          |                                                       |
|---------------------------------------------|----------------------------------------------------------|-------------------------------------------------------|
| <b>Catalog Number:</b><br>84803-2-PBS       | <b>GenBank Accession Number:</b><br>BC026967             | <b>Purification Method:</b><br>Protein A purification |
| <b>Size:</b><br>1 mg/ml                     | <b>GeneID (NCBI):</b><br>25921                           | <b>CloneNo.:</b><br>242192B9                          |
| <b>Source:</b><br>Rabbit                    | <b>UNIPROT ID:</b><br>Q9C0B5                             |                                                       |
| <b>Isotype:</b><br>IgG                      | <b>Full Name:</b><br>zinc finger, DHHC-type containing 5 |                                                       |
| <b>Immunogen Catalog Number:</b><br>AG15022 | <b>Calculated MW:</b><br>715 aa, 78 kDa                  |                                                       |

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01576-2, ZDHHC5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84803-1-PBS. Detection antibody: 84803-2-PBS. Standard: Ag15022. Range: 0.625-80 ng/mL.