

For Research Use Only

SIRT7 Recombinant antibody, PBS Only (Detector)

Catalog Number:84805-2-PBS



Basic Information

Catalog Number: 84805-2-PBS	GenBank Accession Number: NM_016538	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 51547	CloneNo.: 242167A8
Source: Rabbit	UNIPROT ID: Q9NRC8	
Isotype: IgG	Full Name: sirtuin (silent mating type information regulation 2 homolog) 7 (S. cerevisiae)	
Immunogen Catalog Number: AG31794	Calculated MW: 45KD	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01572-1, SIRT7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84805-3-PBS. Detection antibody: 84805-2-PBS. Standard: Ag31794. Range: 3.125-100 ng/mL.



Cytometric bead array standard curve of MP01572-2, SIRT7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84805-1-PBS. Detection antibody: 84805-2-PBS. Standard: Ag31794. Range: 3.125-100 ng/mL.