

For Research Use Only

ACSS3 Recombinant antibody, PBS Only

Catalog Number: 84810-4-PBS



Basic Information

Catalog Number:

84810-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9173

GenBank Accession Number:

BC009317

GeneID (NCBI):

79611

UNIPROT ID:

Q9H6R3

Full Name:

acyl-CoA synthetase short-chain
family member 3

Calculated MW:

686 aa, 75 kDa

Observed MW:

70-75 kDa

Purification Method:

Protein A purification

CloneNo.:

242423C6

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

ACSS3(acyl-CoA synthetase short-chain family member 3, mitochondrial) belongs to the ATP-dependent AMP-binding enzyme family. It activates acetate so that it can be used for lipid synthesis or for energy generation. The deduced 686-amino acid protein contains 4 of 5 motifs characteristic of acyl-CoA synthetases. This protein has 2 isoforms produced by alternative splicing. The full length protein has a transit peptide with 29 amino acids.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 84810-4-RR (ACSS3 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84810-4-PBS in a different storage buffer formulation.



THP-1 cells were subjected to SDS PAGE followed by western blot with 84810-4-RR (ACSS3 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84810-4-PBS in a different storage buffer formulation.