

For Research Use Only

TLR6 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84811-3-PBS



Basic Information

Catalog Number: 84811-3-PBS	GenBank Accession Number: BC111755	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10333	CloneNo.: 241926C4
Source: Rabbit	UNIPROT ID: Q9Y2C9	
Isotype: IgG	Full Name: toll-like receptor 6	
	Calculated MW: 480 aa, 56 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01727-1, Human TLR6 Recombinant Matched Antibody Pair - PBS only. 84811-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard SY02306. 84811-2-PBS was HRP conjugated as the detection antibody. Range: 6.25-400 pg/mL.



The mean TLR6 concentration was determined to be 166.5 pg/mL in Daudi cell extract based on a 1.8 mg/mL extract load.