

For Research Use Only

CYP24A1 Recombinant antibody, PBS Only (Capture)

Catalog Number:84821-1-PBS



Basic Information

Catalog Number:

84821-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17089

GenBank Accession Number:

BC109083

GeneID (NCBI):

1591

UNIPROT ID:

Q07973

Full Name:

cytochrome P450, family 24, subfamily A, polypeptide 1

Calculated MW:

514 aa, 59 kDa

Purification Method:

Protein A purification

CloneNo.:

242142B11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01587-1, CYP24A1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84821-1-PBS. Detection antibody: 84821-2-PBS. Standard: Ag17089. Range: 3.125-100 ng/mL