

For Research Use Only

EZH2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84824-1-PBS



Basic Information

Catalog Number: 84824-1-PBS	GenBank Accession Number: BC010858	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2146	CloneNo.: 242309C4
Source: Rabbit	UNIPROT ID: Q15910	
Isotype: IgG	Full Name: enhancer of zeste homolog 2 (Drosophila)	
Immunogen Catalog Number: AG16789	Calculated MW: 751 aa, 86 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01605-2, EZH2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84824-2-PBS. Detection antibody: 84824-1-PBS. Standard: Ag16789. Range: 0.781-100 ng/mL.