

For Research Use Only

FUT4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84825-1-PBS



Basic Information

Catalog Number:

84825-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC136374

GeneID (NCBI):

2526

UNIPROT ID:

P22083

Full Name:

fucosyltransferase 4 (alpha (1,3)
fucosyltransferase, myeloid-specific)

Calculated MW:

530 aa, 59 kDa

Purification Method:

Protein A purification

CloneNo.:

241888C9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01607-1, FUT4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84825-3-PBS. Detection antibody: 84825-1-PBS. Standard: Eg2194. Range: 0.313-40 ng/mL