

For Research Use Only

CSF3R Recombinant antibody, PBS Only (Capture)

Catalog Number: 84827-1-PBS



Basic Information

Catalog Number:

84827-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_000760.3

GeneID (NCBI):

1441

UNIPROT ID:

Q99062-1

Full Name:

colony stimulating factor 3 receptor
(granulocyte)

Calculated MW:

92 kDa

Purification Method:

Protein A purification

CloneNo.:

242246A10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01609-1, G-CSFR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84827-1-PBS. Detection antibody: 84827-2-PBS. Standard: Eg2821. Range: 0.781-100 ng/mL.