

For Research Use Only

MLF1IP Recombinant antibody, PBS Only (Capture)

Catalog Number: 84847-4-PBS



Basic Information

Catalog Number: 84847-4-PBS	GenBank Accession Number: BC031520	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 79682	CloneNo.: 242194E6
Source: Rabbit	UNIPROT ID: Q71F23	
Isotype: IgG	Full Name: MLF1 interacting protein	
Immunogen Catalog Number: AG3858	Calculated MW: 418 aa, 48 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01624-2, MLF1IP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84847-4-PBS. Detection antibody: 84847-3-PBS. Standard: Ag3858. Range: 0.313-20 ng/mL