

For Research Use Only

# FOXP4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84855-3-PBS



## Basic Information

<b>Catalog Number:</b> 84855-3-PBS	<b>GenBank Accession Number:</b> BC052803	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 116113	<b>CloneNo.:</b> 242356E3
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q8IVH2	
<b>Isotype:</b> IgG	<b>Full Name:</b> forkhead box P4	
<b>Immunogen Catalog Number:</b> AG10047	<b>Calculated MW:</b> 678 aa, 73 kDa	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sandwich ELISA standard curve of MP01621-2, Human FOXP4 Recombinant Matched Antibody Pair - PBS only. 84855-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag10047. 84855-3-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL.



The mean FOXP4 concentration was determined to be 5.4 ng/mL in HepG2 cell extract based on a 1.0 mg/mL extract load.