

For Research Use Only

SLC3A1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84856-1-PBS



Basic Information

Catalog Number:

84856-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9525

GenBank Accession Number:

BC022386

GeneID (NCBI):

6519

UNIPROT ID:

Q07837

Full Name:

solute carrier family 3 (cystine, dibasic and neutral amino acid transporters, activator of cystine, dibasic and neutral amino acid transport), member 1

Calculated MW:

685 aa, 79 kDa

Purification Method:

Protein A purification

CloneNo.:

242370C7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01622-1, SLC3A1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84856-1-PBS. Detection antibody: 84856-3-PBS. Standard: Ag9525. Range: 0.625-80 ng/mL