

For Research Use Only

IL-10RB Recombinant antibody, PBS Only (Detector)

Catalog Number: 84858-1-PBS



Basic Information

Catalog Number:

84858-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC001903

GeneID (NCBI):

3588

UNIPROT ID:

Q08334

Full Name:

interleukin 10 receptor, beta

Calculated MW:

37 kDa

Purification Method:

Protein A purification

CloneNo.:

242302A3

Applications

Tested Applications:

Cytometric bead array, Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Sandwich ELISA standard curve of MP01631-1, Human IL-10RB Recombinant Matched Antibody Pair - PBS only. 84858-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg3446. 84858-1-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL



The mean IL-10RB concentration was determined to be 5.84 ng/mL in THP-1 cell extract based on a 4.40 mg/mL extract load.



Cytometric bead array standard curve of MP01631-1, IL-10RB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84858-2-PBS. Detection antibody: 84858-1-PBS. Standard: Eg3446. Range: 0.156-20 ng/mL