

For Research Use Only

DDX58 Recombinant antibody, PBS Only (Detector)

Catalog Number:84861-2-PBS



Basic Information

| | | |
|---|--|---|
| Catalog Number: 84861-2-PBS | GenBank Accession Number: BC132786 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 23586 | CloneNo.: 242424D4 |
| Source: Rabbit | UNIPROT ID: O95786 | |
| Isotype: IgG | Full Name: DEAD (Asp-Glu-Ala-Asp) box polypeptide 58 | |
| Immunogen Catalog Number: AG18585 | Calculated MW: 925 aa, 106 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01632-2, RIG-1/DDX58 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84861-4-PBS. Detection antibody: 84861-2-PBS. Standard: Ag18585. Range: 0.625-80 ng/mL.