

For Research Use Only

SLC5A5 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84864-1-PBS



Basic Information

Catalog Number: 84864-1-PBS	GenBank Accession Number: BC105047	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6528	CloneNo.: 242516B4
Source: Rabbit	UNIPROT ID: Q92911	
Isotype: IgG	Full Name: solute carrier family 5 (sodium iodide symporter), member 5	
Immunogen Catalog Number: AG19504	Calculated MW: 643 aa, 69 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01665-1, SLC5A5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84864-3-PBS. Detection antibody: 84864-1-PBS. Standard: Ag19504. Range: 0.781-100 ng/mL.