

For Research Use Only

CENPF Recombinant antibody, PBS Only (Capture)

Catalog Number: 84865-3-PBS



Basic Information

Catalog Number: 84865-3-PBS	GenBank Accession Number: NM_016343	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1063	CloneNo.: 242277H7
Source: Rabbit	UNIPROT ID: P49454	
Isotype: IgG	Full Name: centromere protein F, 350/400ka (mitosin)	
Immunogen Catalog Number: AG29665	Calculated MW: 368 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01637-2, CENPF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84865-3-PBS. Detection antibody: 84865-1-PBS. Standard: Ag29665. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01637-1, CENPF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84865-3-PBS. Detection antibody: 84865-2-PBS. Standard: Ag29665. Range: 0.625-80 ng/mL