

For Research Use Only

MAPK10 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84875-2-PBS



Basic Information

Catalog Number: 84875-2-PBS	GenBank Accession Number: BC065516	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5602	CloneNo.: 242364B3
Source: Rabbit	UNIPROT ID: P53779	
Isotype: IgG	Full Name: mitogen-activated protein kinase 10	
Immunogen Catalog Number: AG11785	Calculated MW: 464 aa, 53 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01616-1, JNK Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84875-1-PBS. Detection antibody: 84875-2-PBS. Standard: Ag11785. Range: 0.469-60 ng/mL