

For Research Use Only

CHKA Recombinant antibody, PBS Only (Capture)

Catalog Number: 84881-1-PBS



Basic Information

Catalog Number: 84881-1-PBS	GenBank Accession Number: BC036471	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1119	CloneNo.: 242095B11
Source: Rabbit	UNIPROT ID: P35790	
Isotype: IgG	Full Name: choline kinase alpha	
Immunogen Catalog Number: AG4761	Calculated MW: 439 aa, 52 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01644-2, CHKA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84881-1-PBS. Detection antibody: 84881-2-PBS. Standard: Ag4761. Range: 1.563-200 ng/mL.



Cytometric bead array standard curve of MP01644-1, CHKA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84881-1-PBS. Detection antibody: 84881-3-PBS. Standard: Ag4761. Range: 1.563-200 ng/mL.