

For Research Use Only

# MYCN Recombinant antibody, PBS Only (Capture)

Catalog Number: 84886-1-PBS



## Basic Information

Catalog Number:

84886-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27224

GenBank Accession Number:

BC002712

GeneID (NCBI):

4613

UNIPROT ID:

P04198

Full Name:

v-myc myelocytomatosis viral  
related oncogene, neuroblastoma  
derived (avian)

Calculated MW:

50 kDa

Purification Method:

Protein A purification

CloneNo.:

242406A9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01659-1, MYCN Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84886-1-PBS. Detection antibody: 84886-2-PBS. Standard: Ag27224. Range: 0.781-100 ng/mL.