

For Research Use Only

PLSCR1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84888-2-PBS



Basic Information

Catalog Number: 84888-2-PBS	GenBank Accession Number: BC021100	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5359	CloneNo.: 242372F12
Source: Rabbit	UNIPROT ID: O15162	
Isotype: IgG	Full Name: phospholipid scramblase 1	
Immunogen Catalog Number: AG2172	Calculated MW: 318 aa, 35 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01664-1, PLSCR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84888-2-PBS. Detection antibody: 84888-3-PBS. Standard: Ag2172. Range: 0.781-100 ng/mL.