

For Research Use Only

NUAK1 Recombinant antibody, PBS Only (Detector)

Catalog Number:84893-2-PBS



Basic Information

Catalog Number: 84893-2-PBS	GenBank Accession Number: BC152462	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9891	CloneNo.: 242411A3
Source: Rabbit	UNIPROT ID: O60285	
Isotype: IgG	Full Name: NUAK family, SNF1-like kinase, 1	
Immunogen Catalog Number: AG18832	Calculated MW: 661 aa, 74 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01658-2, ARK5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84893-3-PBS. Detection antibody: 84893-2-PBS. Standard: Ag18832. Range: 0.781-100 ng/mL.