

For Research Use Only

# ARK5 Recombinant antibody, PBS Only

Catalog Number: 84893-5-PBS



## Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
84893-5-PBS	BC152462	Protein A purification
Size:	GeneID (NCBI):	CloneNo.:
1 mg/ml	9891	242411F1
Source:	UNIPROT ID:	
Rabbit	O60285	
Isotype:	Full Name:	
IgG	NUAK family, SNF1-like kinase, 1	
Immunogen Catalog Number:	Calculated MW:	
AG18832	661 aa, 74 kDa	

## Applications

Tested Applications:  
IF/ICC, Indirect ELISA

Species Specificity:  
human

## Background Information

NUAK1 is a member of the human adenosine monophosphate (AMP)-activated protein kinase family that has been identified as a key tumor cell survival factor. It is also named as ARK5, KIAA0537, OMPHK1 and belongs to the protein kinase superfamily. The physiological role of NUAK1 is potentially to suppress glucose uptake through negative regulation of INS signaling in oxidative muscle. This protein has 2 isoforms produced by alternative splicing.

## Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed A431 cells using NUAK1 antibody (84893-5-RR, Clone: 242411F1 ) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 84893-5-PBS in a different storage buffer formulation.