

For Research Use Only

ARK5 Recombinant antibody

Catalog Number: 84893-5-RR



Basic Information

Catalog Number: 84893-5-RR	GenBank Accession Number: BC152462	Purification Method: Protein A purification
Size: 1000 µg/ml	GeneID (NCBI): 9891	CloneNo.: 242411F1
Source: Rabbit	UNIPROT ID: O60285	Recommended Dilutions: IF/ICC 1:200-1:800
Isotype: IgG	Full Name: NUAK family, SNF1-like kinase, 1	
Immunogen Catalog Number: AG18832	Calculated MW: 661 aa, 74 kDa	

Applications

Tested Applications: IF/ICC, ELISA	Positive Controls: IF/ICC : A431 cells,
Species Specificity: human	

Background Information

NUAK1 is a member of the human adenosine monophosphate (AMP)-activated protein kinase family that has been identified as a key tumor cell survival factor. It is also named as ARK5, KIAA0537, OMPHK1 and belongs to the protein kinase superfamily. The physiological role of NUAK1 is potentially to suppress glucose uptake through negative regulation of INS signaling in oxidative muscle. This protein has 2 isoforms produced by alternative splicing.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed A431 cells using NUAK1 antibody (84893-5-RR, Clone: 242411F1) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red).