

For Research Use Only

Cleaved Caspase 3/P17/P19 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84912-1-PBS



Basic Information

Catalog Number:

84912-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC016926

GeneID (NCBI):

836

UNIPROT ID:

P42574

Full Name:

caspase 3, apoptosis-related cysteine
peptidase

Calculated MW:

277 aa, 32 kDa

Purification Method:

Protein A purification

CloneNo.:

241972C6

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00104-4, Human Cleaved Caspase 3/P17/P19 Recombinant Matched Antibody Pair - PBS only. 84912-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard SY02375. 84912-1-PBS was HRP conjugated as the detection antibody. Range: 0.625-40 ng/mL



The mean Cleaved Caspase 3/P17 concentration was determined to be 44.51 ng/mL in LNCaP cell extract based on a 2.00 mg/mL extract load.