

For Research Use Only

TCOF1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84914-2-PBS



Basic Information

Catalog Number: 84914-2-PBS	GenBank Accession Number: BC014559	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6949	CloneNo.: 242201G12
Source: Rabbit	UNIPROT ID: Q13428	
Isotype: IgG	Full Name: Treacher Collins-Franceschetti syndrome 1	
Immunogen Catalog Number: AG1447	Calculated MW: 97 kDa, 144 kDa, 152 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01672-2, TCOF1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84914-2-PBS. Detection antibody: 84914-1-PBS. Standard: Ag1447. Range: 0.625-80 ng/mL