

For Research Use Only

GPR56 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84915-3-PBS



Basic Information

Catalog Number:

84915-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001145771.2

GeneID (NCBI):

9289

UNIPROT ID:

Q9Y653

Full Name:

G protein-coupled receptor 56

Calculated MW:

78 kDa

Purification Method:

Protein A purification

CloneNo.:

242249G6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01673-2, GPR56 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84915-1-PBS. Detection antibody: 84915-3-PBS. Standard: Eg2478. Range: 0.781-50 ng/mL.